

Easy
IFU

EasyIFU

User Manual

Version: 1.0

Date: May 2025

Author: Easy Medical Device

Table des matières

Table des matières	2
Introduction.....	3
1. Getting Started	4
2. Managing Your Profile and Company.....	4
3. Creating and Managing Projects	7
4. eIFU Management Workflow	8
5. Label Management Workflow	16
6. Production Management	22
7. Notifications and Communication.....	25
8. Logs and History	25
9. Role-Based Access Overview	26
Conclusion	26



Introduction

Welcome to the EasyIFU User Manual! This guide is designed to help users effectively navigate and use the EasyIFU platform for managing electronic Instructions for Use (eIFUs), product labels, and production tasks.

What is EasyIFU?

EasyIFU is a cloud-based solution tailored for medical device companies to:

- Create and manage eIFUs and regulatory labels
- Ensure proper review, approval, and release workflows
- Automatically generate multilingual QR codes for public access
- Facilitate dynamic label printing and audit tracking

With user roles like Creator, Approver, Releaser, Producer, and OrgAdmin, EasyIFU offers controlled access and traceability across the entire documentation lifecycle.

1. Getting Started

Once you are logged in, you will be directed to the Dashboard, which provides an overview of your projects, status updates, and shortcuts to key modules.

The dashboard provides a comprehensive overview of the organization's operations:

- Welcome back, Org Admin!** Here's your dashboard overview. Villa N°3 Rue De La Guinee Dior Jamaa Rabat RABAT Morocco
- Label Status Over Time:** A line chart showing the number of labels over time, categorized by status: Labels Created (green), Labels Approved (blue), Labels Released (orange), and Labels Rejected (red).
- Key Metrics:**
 - Ongoing Labels: 0
 - Received Labels: 0
 - Total Labels: 0
 - Total Orders Printed: 0
- Received Labels:** A table showing received labels with columns: ID, Order Name, Status, and Open. Status is currently Open.
- Received eIFUs:** A table showing received eIFUs with columns: ID, Order Name, Status, and Open. Status is currently Open.
- Label Status Distribution:** A bar chart showing the distribution of labels by status: rejected, released, pending_release, Produced, approved, pending_approval, and Draft.
- eIFU Count by Status:** A bar chart showing the count of eIFUs by status: rejected, released, pending_release, Produced, approved, pending_approval, and Draft.
- Monthly Metrics Comparison:** A chart comparing monthly metrics across different categories.
- eIFU Trends Over Time:** A line chart showing trends over time for eIFU creation.
- Recent Labels:** A table showing recent labels with columns: ID, Name, Version, Status, and Open.

2. Managing Your Profile and Company

- Go to your **profile** icon in the top-right corner.

Welcome back, Org Admin! Here's your dashboard overview.
Villa N°3 Rue De La Guinee Dior Jamaa Rabat RABAT Morocco

0 Ongoing Labels	0 Received Labels
0 Total Labels	0 Total Orders Printed

Label Status Over Time:

Received Labels: [Show More ▶](#)

ID	Order Name	Status	Open
There are currently no labels pending for approval or release.			

Received eIFU's: [Show More ▶](#)

ID	Order Name	Status	Open
No Instructions for Use (IFUs) are awaiting approval or publication at this time.			

- Click on "Profile" to update your personal information.

-/ Dashboard / Account

First name: Org
Last name: Admin
Email: loubnaelghali889@gmail.com
Role: Admin, Creator, Release, Approver, Producer,
Status: active
Creation Date: 19/05/2025

Hide Form

First Name: Org Last Name: Admin
Password
User Roles: Admin Creator Approver Release Producer

Update

- Navigate to "My Company" to view your company details.

Welcome back, Org Admin! Here's your dashboard overview.
Villa N°3 Rue De La Guinee Dior Jamaa Rabat RABAT Morocco

0 Ongoing Labels	0 Received Labels
0 Total Labels	0 Total Orders Printed

Label Status Over Time:

Received Labels: [Show More ▶](#)

ID	Order Name	Status	Open
There are currently no labels pending for approval or release.			

Received eIFU's: [Show More ▶](#)

ID	Order Name	Status	Open
No Instructions for Use (IFUs) are awaiting approval or publication at this time.			

- Click on "**Create User**" to create new users by entering their name, email, and password, and assigning the appropriate role(s).

- Click on "**Users**" to manage users: activate/deactivate/ delete, and update roles.

#	Full Name	Role	Created	Email	Status	Manage
1	Org Admin	Admin, Creator, Release, Approver, Producer,	19/05/2025	loubnaelghali889@gmail.com	Active	
2	Approver Releaser	Approver, Release,	19/05/2025	lubnaetude@gmail.com	Active	

3. Creating and Managing Projects

Follow these steps to create and manage projects:

- Ensure your user role is set to Creator. Only users with this role can create projects.
- Go to the left-hand menu and click on the "Projects" tab.

- Click the "New Project" button.

- Fill in the Project Name and Description in the fields provided.

No Projects Created Yet

Project Information

1- Project Name:

2- Project Description:

Start Creating **Close**

- Click "Start Creating" to finalize the new project setup.

Project created successfully

#	Project Name	Description	CreatedAt	Details	Delete
1	Easy Project	Lorem Ipsum has been...	19/05/2025 12:24:00		

Note: You can only delete a project if it has no products or labels associated with it.

4. eIFU Management Workflow

- Go to the eIFU section.

#	Project Name	Description	CreatedAt	Details	Delete
1	Easy Project	Lorem Ipsum has been...	19/05/2025 12:24:00		

The screenshot shows the Easy IFU application interface. On the left is a sidebar with the 'Org Admin' user information and navigation links for Home, Projects, Orders, and contact Us. The main content area has a header with 'Products' and 'Electronic Instructions for Use'. Below the header is a breadcrumb path '- / Dashboard / Products'. A red arrow points to the 'Electronic Instructions for Use' tab. The main content area contains a table with columns: #, Product Name, Description, CreatedAt, and Details. A message at the bottom says 'No product created'.

- Create an eIFU product by clicking on "New eIFU" and fill in the Name in the field provided.

The screenshot shows the Easy IFU application interface. On the left is a sidebar with the 'Org Admin' user information and navigation links for Home, Projects, Orders, and contact Us. The main content area has a header with 'Products' and 'Instructions for Use'. Below the header is a breadcrumb path '- / Dashboard / Instructions For Use'. A red arrow points to the 'New eIFU' button. The main content area contains a table with columns: #, Instruction Name, CreatedAt, and Details. A message at the bottom says 'No Electronic Instructions for Use created yet'.

The screenshot shows a modal dialog box titled 'New eIFU'. Inside the dialog, there is an input field labeled 'Instructions Name:' containing the value 'Easy eIFU'. At the bottom right of the dialog are two buttons: 'Create' and 'Close'.

- For your first eIFU, you should create a prefix that will be used in the eIFU's URL. This ensures all your published instructions follow a consistent naming pattern.

The screenshot shows the Easy IFU dashboard with a modal dialog centered over the main content. The modal has a white background with a large blue circular icon containing a white letter 'i'. Below the icon, the text "Instruction URL" is displayed. A message below states: "Instruction URL prefix is missing. Please create one to proceed." At the bottom of the modal are two buttons: a blue "Create" button and a red "Cancel" button.

#	Instruction Name	CreatedAt	Details
1	Easy eIFU	20/05/2025 15:55:12	

The screenshot shows a configuration page titled "Set Your eIFU Prefix". It features a text input field with the placeholder "Choose a prefix for your instruction URL" and the value "EASY" entered. Below the input field is a "Save Prefix" button. A link "Want to do this later? Go to Dashboard" is visible at the bottom of the form. An example URL "Example: easyifu.com/EASY-1234" is shown at the bottom right.

- Now, you can create an eIFU.

The screenshot shows the dashboard after creating a new eIFU. The table now contains one row with the following data:

#	Instruction Name	CreatedAt	Details
1	Easy eIFU	20/05/2025 15:55:12	

The screenshot shows the 'Org Admin' section of the software. On the left, there's a sidebar with links to Home, Projects, Orders, and contact Us. A yellow 'Upgrade to Premium' button is at the bottom. The main area shows a table for 'elfUs' with columns: #, Created At, elfU Name, Description, Versions, Status, langue, Edit, and Details. A red arrow points to the 'New elfU' button at the top left of the table area.

- Fill in the elfU name and description, upload the IFU document, select the language, and complete the checklist. Then save the instruction.

This screenshot shows the 'Instructions for Use Form' creation page. On the left, there's a preview of the document content, which includes regulatory references (ISO 13485:2016, ISO 20417:2021), a user manual section, and a general provisions section. On the right, there are fields for 'Name' (Easy elfU), 'Description' (a placeholder text about the document providing safe instructions), 'Select Language' (Global), and a 'Checklist' section with one question. A red arrow points to the 'File Selected' field.

This screenshot shows the 'Instructions for Use Form' creation page after filling out the checklist. The checklist contains several questions with radio button options. At the bottom right, there is a large blue 'save instruction' button with a red arrow pointing to it.

- The eIFU has been created. You need to go back for the review.

No PDF selected

Drag & drop your PDF here or click to upload

Name: _____

Description: _____

Select Language
Language selected: _____

Checklist:

- Now, it's time to review and publish the eIFU.

New eIFU

#	Created At	eIFU Name	Description	Versions	Status	langue	Edit	Details
1	20/05/2025 22:18:49	Easy eIFU	This document provid...	0.1	Draft	Global		

Version: 1

Regulatory References

ISO 13485:2016	4.2.3
ISO 20417:2021	All

User Manual [enter Product Name]

Version: [enter content] Date: [enter content]
[Placeholder Company Address] Email: [enter content] Website: [enter content]
[Placeholder CE Sign] [Placeholder UDI]

[enter Product Name] is a class [enter content] medical device in accordance with 93/42/EEC (for MDD) // Regulation (EU) 2017/745 (for MDR). UMDNS Code: [enter content]

Symbols:
Safety Information: [Placeholder Symbol] This symbol indicates important safety-related information.
Additional Information: [Placeholder Symbol] This symbol indicates supplementary information.

On your introductory pages, explain the symbols that you use in the following sections. Take into account standards like ISO 15223.

General Requirements: Instructions for use are regulated in Section 23 of the General Requirements (MDD) and Chapter 3, Section 23 of the General Safety and Performance Requirements (MDR). You may be able to make a case for why your users must not be provided with instructions for use, based on the exception stated in Section 13.1 (MDD) or Section 23.1.d (MDR):

"By way of exception, no such instructions for use are needed for devices in Class I or IIa if they can be used safely without any such instructions."

This document provides detailed instructions for the safe and effective use of the Easy eIFU system. It includes step-by-step guidance, safety information, and compliance requirements to support users throughout the product lifecycle.

Questions

True	False	Ap
the name, registered trade name or registered trade mark of the manufacturer and the address of its registered place of business ?	✓	
Indication if your device contain a medicinal substance, including a human blood or plasma derivative ?	✓	
Warnings or precautions to be taken ?	✓	
if the device is intended for single use, an indication of that fact ?		
Indication if your device contains CMR: List of substances ?		
If device sterile, an indication of its sterile state and the sterilisation method ?	✓	
Indication for any special storage and/or handling conditions ?	✓	

choose user

- Choose and send the eIFU to the approver (Note: the creator of the eIFU cannot be an approver).

The screenshot shows a web-based application for managing medical device information. At the top, there are navigation icons and a status bar indicating the document ID (4a81b920-61a...), page number (1 / 4), zoom level (90%), and other document details. A modal dialog box titled "choose user" is open, listing two users: "Org Admin" and "Approver Releaser". Both users have checkboxes next to them labeled "Status: Draft" and "Created by: Admin Org". Below the list is a "Close" button. To the right of the modal, the main content area displays a "User Manual [enter Product Name]" section, regulatory references (ISO 13485:2016, ISO 20417:2021), and various product details. On the far right, there are "True", "False", and "Ap" columns for different questions.

This screenshot shows the same interface after the user has selected the approver. The "Approver Releaser" is now listed under "Received" at the bottom of the screen. A large blue "Send" button is prominently displayed in the center. A red arrow points to this button, indicating the next step. The rest of the interface remains largely the same, showing the user manual, regulatory references, and product details.

This final screenshot shows the completed eIFU document. The "Approver Releaser" is listed under "Received". The "Send" button is no longer present. Instead, a red arrow points to a blue "←" button located in the top left corner of the main content area. The "Description" section is visible, containing a brief overview of the document's purpose. The "Questions" section is also shown, with the "True", "False", and "Ap" columns for the previously listed questions.

Org Admin
loubnaelghali889@gmail.com

New eIFU Global

#	Created At	eIFU Name	Description	Versions	Status	langue	Edit	Details
1	21/05/2025 00:03:44	Easy eIFU	This document provid...	0.1	pending to approve	Global		

Home Projects Orders contact Us

Upgrade to Premium

Version: 1

- For the Approver / Releaser:

Approver Releaser
lubnaelghali889@gmail.com

Welcome back, Approver Releaser! Here's your dashboard overview.
Villa N°3 Rue De La Guinee Dior Jamaa Rabat RABAT Morocco

Ongoing Labels: 0 Received Labels: 0

Total Labels: 0 Total Orders Printed: 0

Label Status Over Time:

Received Labels:
Show More ►

ID	Order Name	Status	Open
ID	Order Name	Status	Open

Approver Releaser
lubnaelghali889@gmail.com

Received Labels Received eIFU's

Received Projects: 0

#	Label Name	Description	Manage
---	------------	-------------	--------

You have no Label to approve or release yet

Home Projects contact Us

Upgrade to Premium

Version: 1

Approver Releaser
lubnaetude@gmail.com

Received Labels **Received IFUs**

Received Project: 1

#	IFU Name	Description	Manage
1	Easy eIFU	This document provid...	

Home **Projects** **Contact Us**

Upgrade to Premium

Version: 1

Easy eIFU

You may accept or reject it

Accept **Reject**

Description :

This document provides detailed instructions for the safe and effective use of the Easy eIFU system. It includes step-by-step guidance, safety information, and compliance requirements to support users throughout the product lifecycle.

Language: Global
 Version: 0.1
 Status: pending to approve
 Created by: Admin Org

Questions

the name, registered trade name or registered trade mark of the manufacturer and the address of its registered place of business ?

Indication if your device contain a medicinal substance, including a human blood or plasma derivative ?

Easy eIFU

Send the eIFU to the Publisher.

-- select releaser -- **Send to Publisher**

Description :

This document provides detailed instructions for the safe and effective use of the Easy eIFU system. It includes step-by-step guidance, safety information, and compliance requirements to support users throughout the product lifecycle.

Language: Global
 Version: 0.1
 Status: approved
 Created by: Admin Org
 Approved by: Releaser Approver

Questions

the name, registered trade name or registered trade mark of the manufacturer and the address of its registered place of business ?

Approver Releaser
lubaetaude@gmail.com

Home Projects contact Us

Upgrade to Premium

Version: 1

- The eIFU is published

Regulatory References

ISO 13485:2016	4.2.3
ISO 20417:2021	All

User Manual [enter Product Name]

Version: [enter content] Date: [enter content]
[Placeholder Company Address] Email: [enter content] Website: [enter content]
[Placeholder CE Sign] [Placeholder UDI]
[enter Product Name] is a class [enter content] medical device in accordance with 93/42/EEC (for MDD) // Regulation (EU) 2017/745 (for MDR). UMDNS Code: [enter content]

Symbols:
Safety Information: [Placeholder Symbol] This symbol indicates important safety-related information.
Additional Information: [Placeholder Symbol] This symbol indicates supplementary information.

On your introductory pages, explain the symbols that you use in the following sections. Take into account standards like ISO 15223.

General Provisions: Instructions for use are regulated in Section

Description :

This document provides detailed instructions for the safe and effective use of the Easy eIFU system. It includes step-by-step guidance, safety information, and compliance requirements to support users throughout the product lifecycle.

Language: Global
Version: 1.0
Status: published
Created by: Admin Org
Approved by: Releaser Approver
Published by: Releaser Approver

Download QR Code

QR Code: easyifu.com/EASY17

Questions

the name, registered trade name or registered trade mark of the manufacturer and the address of its registered place of business ? ✓

Indication if your device contain a medicinal substance, including a human blood or plasma derivative ?

Indication for any special storage and/or handling conditions ?

5. Label Management Workflow

Create Label Product:

- Answer 11 sets of regulatory questions (Manufacturer, Legislation, Product Info, etc.)
- Provide Intended Purpose in multiple languages
- Define IFU link or QR code source

Product Information

1- Product Name:
Easy Product Test

2- Product Description:
Lorem Ipsum has been the industry's standard dummy text ever since the 1500s, w...

Start Creating **Close**

#	Product Name	Description	CreatedAt	Details
1	Easy Product Test	Lorem Ipsum has been...	20/05/2025 15:54:27	Continue

#	Product Name	Description	CreatedAt	Details
1	Easy Product Test	Lorem Ipsum has been...	20/05/2025 15:54:27	Continue

Manufacturer Information

1- Do You have a Distributor ?
Yes: No:

2- European Authorized Representative information:

European Authorized Representative Name:

European Authorized Representative Address:

Importer Name:

Importer Address:

- Answer the questions in each step according to your product. Step 7 is where you should add your eIFU. If you already have one in EasyIFU, you can select it; otherwise, leave it empty and add it later.

The screenshot shows the Easy IFU application interface. On the left, there's a sidebar with 'Org Admin' and 'loubnaelghali889@gmail.com'. Below that are links for 'Home', 'Projects', 'Orders', and 'contact Us', along with a yellow 'Upgrade to Premium' button. The main area has a header 'All Labels Draft Approved Released Rejected' and a sub-header '~/Dashboard / Labels'. A red box highlights the 'New Label' button. Below it is a table with columns: #, Label Name, Version, Description, Status, Details, Edit, and Created At. The message 'No Projects Created' is displayed.

Generate Label:

- Select a template
- Preview and edit the label

This screenshot shows a modal dialog titled 'Project Information'. Inside, there are three input fields: '1- Label Name:' containing 'Easy Label Lot', '2- Label Description:' containing 'UltraSafe is a single-use, retractable hypodermic needle designed to minimize the', and '3- choose the intended purpose for this label:' with checkboxes for 'English' and 'Swedish' both checked. At the bottom are two buttons: 'Start Creating' on the left and 'Close' on the right.

This screenshot shows the 'Labels' table from the previous screen. A new row has been added with the following data: #1, Label Name 'Easy Label Lot', Version '0.1', Description 'UltraSafe is a singl...', Status 'Draft', and Created At '21/05/2025 00:37:17'. A red arrow points to the 'Edit' button in the 'Edit' column for this row.

#	Label Name	Version	Description	Status	Details	Edit	Created At
1	Easy Label Lot	0.1	UltraSafe is a singl...	Draft			21/05/2025 00:37:17

- To choose between templates

Easy Label Lot

size: 150mmx100mm

Description :

UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Version: 0.1
Created by:OrgAdmin

Show Logs

- To update the label

Easy Label Lot

size: 150mmx100mm

Description :

UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Version: 0.1
Created by:OrgAdmin

Show Logs

- Send it to the approver

Easy Label Lot

size: 150mmx100mm

Description :

UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Version: 0.1
Created by:OrgAdmin

Show Logs

- Approval Workflow :

Received Labels Received IFUs

Received Projects: 0

#	Label Name	Description	Manage
1	Easy Label Lot	UltraSafe is a singl...	

Approver Releaser lubnaetude@gmail.com

Home Projects contact Us

Upgrade to Premium Version: 1

Activate Windows Go to Settings to activate Windows.

Approver Releaser lubnaetude@gmail.com

UltraSafe Hypodermic Needle

Size: 150mmx100mmmm

UltraSafe is a non-invasive continuous glucose monitoring (CGM) system designed to track real-time blood glucose levels in diabetic patients. It helps manage diabetes by providing trend analysis and alerts for low blood sugar levels. The device is intended for use in conjunction with medical health applications.

The UltraSafe Hypodermic Collection Needle is designed for safe and effective blood sampling and collection in healthcare settings, including hospitals, clinics, and laboratories. It is intended for single-use only and ensures a smooth and reliable blood draw while minimizing patient discomfort.

XZ Medical Solutions
Villa N°1 Rue De La Gelinie Dier
RABAT
Morocco

EU Medical Solutions
EU REP 458 Regulatory Rd, Floor 5,
Brussels, Belgium, B-1000

Marketing Services Ltd.
ET Regulatory Lab, Madrid,
Spain

2027-04-21 2027-mm-dd

REF USH-1001 LOT XXXXXX

CE 2797 UK CR0086

UDI

You may accept or reject it.

Accept Reject

Description :
UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Version: 0.1 Created by: Admin Org

Activate Windows Go to Settings to activate Windows.

- If rejected, edit and resubmit

All Labels Draft Approved Released Rejected

~/Dashboard / Labels

New Label

#	Label Name	Version	Description	Status	Details	Edit	Created At
1	Easy Label Lot	0.1	UltraSafe is a singl...	rejected			21/05/2025 00:37:17

Org Admin lubnaelghali809@gmail.com

Home Projects Orders contact Us

Upgrade to Premium Version: 1

Activate Windows Go to Settings to activate Windows.

Org Admin loubnaelghali889@gmail.com

Home Projects Orders contact Us Upgrade to Premium Version: 1

All Labels Draft Approved Released Rejected

/ Dashboard / Labels

New Label

#	Label Name
1	Easy Label Lot

Please be aware that updating certain information will require creating a new product with a new Unique Device Identification (UDI) code. Changes have been made to the following aspects :

- Name or trade name of the device
- Device version or model
- Labelled as single use
- Packaged sterile
- Need for sterilization before use
- Quantity of devices provided in a package
- Critical warning or contra-indications

Continue **Close**

Org Admin loubnaelghali889@gmail.com

Update Label Templates

Back Print

size: 150mmx100mm

Easy Label Lot

Send Project To: Choose User: Send ➔

UltraSafe Hypodermic Needle

Description :

UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Version: 0.2 Created by: OrgAdmin

Show Logs

Activate Windows Go to Settings to activate Windows.

UltraSafe Hypodermic Needle

XYZ Medical Solutions Villa De La Guine Dior
Jalil Ben Rabat
RABAT
Morocco

Medley Services Ltd.
67 Regentary Lane, Madrid,
Spain

2027-04-21

7777-mn-dd

LOT XXXXXX

REF USN-1001

UltraSafe-22G

EU MedReg Solutions
REF 456 Registry Rd, Floor 5,
Brussels, Belgium, B-1000

CE 2797

UK REP

UDI

easylabel.com/EASYIFU

21/05/2023

- Once approved, send to Releaser
- Releaser publishes label

Org Admin loubnaelghali889@gmail.com

Home Projects Orders contact Us Upgrade to Premium Version: 1

All Labels Draft Approved Released Rejected

/ Dashboard / Labels

New Label

#	Label Name	Version	Description	Status	Details	Edit	Created At
1	Easy Label Lot	1.0	UltraSafe is a singl...	released			21/05/2025 00:50:24

Activate Windows Go to Settings to activate Windows.

IFU Description Creation - Google Chrome

6. Production Management

Users with the **Producer** role can:

- Create a new order
- Select released labels
- Enter dynamic data (e.g., lot number, serial number)
- View and print the final label
- Confirm the printing status

Welcome back, Org Admin! Here's your dashboard overview.
Villa N°3 Rue De La Guinee Dior Jamaa Rabat RABAT Morocco

Ongoing Labels: 0

Received Labels: 0

Total Labels: 2

Total Orders Printed: 0

Label Status Over Time:

Number of labels
0.9
0.8
0.7
0.6
0.5
0.4
0.3
0.2
0.1
0
5/2025

Labels Created (blue bar)
Labels Approved (light blue bar)
Labels Released (yellow bar)
Labels Rejected (pink bar)

Received Labels: [Show More ▶](#)

ID	Order Name	Status	Open
There are currently no labels pending for approval or release.			

Received eIFU's: [Show More ▶](#)

ID	Order Name	Status	Open
No Instructions for Use (IFUs) are awaiting approval or publication at this time.			

Label Status Distribution:

eIFU Count by Status:

Activate Windows
Go to Settings to activate Windows.

Version: 1
<https://www.easyifu.com/dashboard>

-/ Dashboard / Orders

New Order

No Order Created

LOT number Serial number [search by LOT number](#)

Activate Windows
Go to Settings to activate Windows.

Version: 1

Dashboard / Templates

#	Template Name	Description	createdAt	Details
1	Easy Label Lot	UltraSafe is a single-use...	21/05/2025 00:50:24	

search by name

Home Projects Orders contact Us

Upgrade to Premium

Version: 1

Activate Windows
Go to Settings to activate Windows.

Back

Easy Label Lot

UltraSafe Hypodermic Needle

Dynamic Data:

Prefix:

Suffix: 0001

Lot Number: Prefix + Suffix

Manufacturer Date: 2025-05-21

how many label you want to print? : 1

Save Order

Label Name: Easy Label Lot

Label description: UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to biohazardous pathogens.

Activate Windows
Go to Settings to activate Windows.

Back

size: 150mmx100mm


UltraSafe Hypodermic Needle

UltraSafe Hypodermic Needle


UltraSafe Hypodermic Needle

Activate Windows
Go to Settings to activate Windows.

- Setup your printer and launch the print



- Then confirm if everything was good



[← Back](#)

Date: 21/05/2025, 12:56:06 AM
Printed by: Org Admin
3 label(s) Printed

UltraSafe Hypodermic Needle MD

SV smartGluco is a new invasive continuous glucose monitoring (CGM) system designed to track and store blood glucose levels in diabetic patients. It helps manage diabetes by providing trend analysis, predictive alerts for hypo- and hyperglycemia, and seamless integration with mobile health applications.

BD UltraSafe Venous Blood Collection Needle is designed for safe and efficient venous blood collection in healthcare settings, including hospitals, clinics, and laboratories. It is intended for single-use only and ensures a smooth and reliable blood draw while minimizing patient discomfort.

VZV Medical Solutions Villa N° 3 Rue De La Guinée Hor Jamaa Rabat Morocco	ModReg Services Ltd. 67 Regulatory Lane, Madrid, Spain
	
2027-04-21	2025-05-21
REF USN: 1001	LOT EMD1
	
EU REF: 450 Regulatory Rd, Floor 5, Brussels, Belgium, B-1000	UK REF:
	
cepsys.com/000117	CE 2797
UDI	UK CR 0086
	101123456789012345671112345672317123456231234565 LAL: 00000001 Y1.0 21/05/2025

Print Info :

Label Size: 150mmx100mm
 Label Id: LBL-0000001-V1.0

Other Info :

Label Name: Easy Label Lot
 Version: 1.0
 Label description: UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Show Logs ↗

Activate Windows
Go to Settings to activate Windows.

7. Notifications and Communication

- Automatic email notifications are sent when:
 - A user account is created or activated
 - A task (approval/release) requires your action
 - Payments or changes in company status occur
 - Use the "**Contact Us**" form for inquiries
 - Submit feedback via the "**Feedback**" tab

8. Logs and History

- Access logs for all actions (e.g., project edits, label submissions)
 - Download audit trails for compliance

[← Back](#)

Date: 21/05/2025, 12:56:06 AM
Printed by: Org Admin
3 label(s) Printed

UltraSafe Hypodermic Needle MD

SV:SmartGluco is a non-invasive continuous glucose monitoring (CGM) system designed to track and store blood glucose levels in diabetic patients. It helps manage diabetes by providing trend analysis, predictive alerts for hypo- and hyperglycemia, and seamless integration with mobile health applications.

UltraSafe Hypodermic Blood Collection Needle is designed for safe and efficient venous blood collection in healthcare settings, including hospitals, clinics, and laboratories. It is intended for single-use only and ensures a smooth and reliable blood draw while minimizing patient discomfort.

YVZ Medical Solutions
Villa N° 3 Rue De La Guinée Hor
Jaume Rabat
IBARAK
Morocco

MedReg Services Ltd.
67 Regent Lane, Madrid,
Spain

2025-05-21

2027-04-21

LOT EMD1

REF USH-1001

REF UltraSafe-226

EU MedReg Solutions
49A Regentbury Rd, Floor 5,
Brussels, Belgium, B-1000

UK REP

UDI

CE 2797

UKCA 0086

8011234567890123451112345622117721984231818805
LBL_00000001_V1.0 21/05/2025

Print Info :

Label Size: 150mmx100mm
 Label Id: LBL-00000001_V1.0

Other Info :

Label Name: Easy Label Lot
 Version: 1.0
 Label description: UltraSafe is a single-use, retractable hypodermic needle designed to minimize the risk of needlestick injuries and enhance injection safety in medical settings. The device incorporates an automatic safety mechanism that retracts the needle into the syringe barrel after use, preventing accidental exposure to bloodborne pathogens.

Show Logs ↗

Activate Windows
Go to Settings to activate Windows.

9. Role-Based Access Overview

Role	Key Permissions
Creator	Create/edit projects, labels, and eIFUs
Approver	Approve/reject eIFUs and labels
Releaser	Publish eIFUs and labels
Producer	Finalize and print labels for production
Admin	Manage users and company settings

Conclusion

You are now ready to efficiently manage your EasyIFU operations. For technical support or training, contact the EasyIFU support team directly through the Help section.